



Hillrom™

Welch Allyn® RetinaVue® 700 Pedestal

Convert the handheld RetinaVue 700 Imager into a tabletop format with ease using the RetinaVue 700 Pedestal.

- Multiple adjustment points for flexible patient positioning
- Quick release enables “take it and go” portability for mobile exams
- Aids in stabilizing across all patient populations
- Minimal to no training required across staff
- Ideal for high-volume clinics



PRODUCT SPECIFICATIONS

| | |
|---------------|-----------------|
| Weight | 8.4lbs (3.8kg) |
| Height | 5in (12.8cm) |
| Length | 14.5in (36.8cm) |
| Width | 11in (27.9cm) |
| Raising Angle | 70° |

ORDERING INFORMATION

| | |
|--------|--|
| 728134 | Welch Allyn Pedestal Base for RetinaVue 700 Imager |
|--------|--|



For more information, please visit retinavue.com.

hillrom.com

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2020 Welch Allyn, Inc. ALL RIGHTS RESERVED. APR84702 Rev 1 80023924 Ver A 20-MAY-2020 ENG – US